

# Human Adrenal gland tumor Cell lines

Part of the CLS cell bank

CLS Cell Lines Service



**Table 1: Human Adrenal gland cancer cell lines: Origin and General Characteristics**

Name of cell line	Cell type	Organism, Ethnicity	Age / Gender	Tissue / Disease	Morphology	Growth properties	CLS order no.
NCI-H295R <sup>1</sup>	Human Adrenal gland carcinoma cell line	Homo sapiens (Human) / Caucasian	48 / Female	Adrenal gland / carcinoma	Epithelial	Monolayer, adherent	300483
SW-13 <sup>2</sup>	Human Adrenal gland carcinoma cell line	Homo sapiens (Human) / Caucasian	55 / Female	Adrenal gland / small cell carcinoma	Epithelial	Monolayer, adherent	300349

**Table 2: Human Adrenal gland cancer cell lines: Special Features**

Name of cell line	Cell type	Cell Marker	Virus Susceptibility	Mutations	Secretion of Products	Ref ID in Cellosaurus <sup>3</sup>	CLS order no.
NCI-H295R <sup>1</sup>	Human Adrenal gland carcinoma cell line				Aldosterone; cortisol; C19 steroids	RRID:CVLL_0458	300483
SW-13 <sup>2</sup>	Human Adrenal gland carcinoma cell line	Isoenzymes: G6PD, B	Vesicular stomatitis (Indiana); poliovirus 1		Reverse Transcriptase:neg	RRID:CVCL_0542	300349

Information on cell culture conditions, authentication data and others can be found on the website: [www.clsgmbh.de](http://www.clsgmbh.de)

## References:

1. Gazdar AF et al. Establishment and characterization of a human adrenocortical carcinoma cell line that expresses multiple pathways of steroid biosynthesis. Cancer Res 50: 5488-96, 1990.
2. Lasfargues EY et al. Cultivation of human breast carcinomas. J Natl Cancer Inst 21: 1131-47, 1958.
3. <http://web.expasy.org/cellosaurus/> - the cellosaurus represents a detailed data collection bank of a plethora of cell line relevant data from various cell banks.

All of the products listed in Table 1/Table 2 are intended for research use only, not for use in human, therapeutic or diagnostic applications.

The General Terms and Conditions of Supply of CLS Cell Lines Service GmbH are valid. According to the Terms, the products are not intended to be resold and/or modified for resale. Their usage in order to provide commercial services or to manufacture commercial products is prohibited unless approved in writing by CLS.

Please contact service (a) [clsgmbh.de](mailto:info@clsgmbh.de) if you have further questions or concerns.