

Karyotype: SW-872

Cell type:	Human Liposarcoma cell line
CLS order numbers:	Cryovial: 300405, Proliferating culture: 330405
Organism / Ethnicity:	<i>Homo sapiens</i> (human) / Caucasian
Age:	36 years
Gender:	Male
Morphology:	Fibroblast
Growth Properties:	Adherent
Passage:	37
Lot #:	300405-710
Banding technique:	G-banding
Metaphase cells karyotyped	5
Short ISCN:	75,XY, (range = 69-79), N1, N2, N10, N18, N22 = 2 Copies N15, N16, N19, N21 = 3 Copies, N3, N4, N5, N6, N7, N8, N9, N11, N12, N13, N14, N17 = 4 Copies, N20 = 5 Copies

